

Waters™



Pipette+

Qualification Procedure



Contents

About this manual	3
1. Installation Qualification (IQ).....	3
1.1. IQ procedure	3
A) User manuals and documents.....	4
B) Installation of Pipette+.....	4
C) Power On the device	4
D) Placement of pipettes on Pipette+ rack.....	5
E) Pairing of Pipette+ to OneLab	5
1.2. IQ test report.....	7
2. Operational Qualification (OQ).....	8
2.1. OQ test plan	8
2.2. Getting started	9
A) Experiment execution and verification correct execution	9
B) Collection of Logfiles	12
2.3. OQ test report	12
3. Contact	13

About this manual

The procedure defines the general steps that should be performed to ensure that Pipette+ is installed and functioning properly.

1. Installation Qualification (IQ)

Document #: 618.002		Revision: 7	
	Function	Name	Date
Author	Lab Automation Service & Support Engineer	L.M. Turrado	25 Feb 2022
Reviewed by	Global Service & Support Sr Manager	N. Vamvoukas	25 Feb 2022

The Installation Qualification (IQ) is documented proof that equipment has been delivered and installed in accordance with the requirements and statutory safety regulations stipulated in the design qualification. The documentation for an Installation Qualification consists of:

- IQ test plan
- IQ report

1.1. IQ procedure

The **IQ procedure** contains detailed information of the tests to be performed on the respective equipment in the context of the Installation Qualification. The IQ procedure is composed of the following series of basic checks:

- Verification that all documents required for the correct use of Pipette+ and its software are given to the customer.
- Verification that the environmental conditions where Pipette+ is installed fulfills the requirements present in the “Quick Start Guide”.
- Verification of the correct installation of Pipette+.

After completion of the Installation Qualification, the results are summarized, evaluated and compiled in an IQ report.

The Installation Qualification template for Pipette+ is reported in the next page.



INSTALLATION QUALIFICATION			
System name:		Pipette+	
Serial number:		618. _____	
Document #	618.002	Revision	7

A) User manuals and documents

Verification that all the documentations and manuals are correctly delivered		
Pipette+ Quick Start Guide	OK <input type="checkbox"/> NOK <input type="checkbox"/>	Comment (if none, check this checkbox <input type="checkbox"/>)

B) Installation of Pipette+

Pipette+ Installation		
Verify that the bench on which Pipette+ is placed is stable, flat enough and dry.	OK <input type="checkbox"/> NOK <input type="checkbox"/>	Comment (if none, check this checkbox <input type="checkbox"/>)

C) Power On the device

Installation of the power line.		
Verify that the power supply provided by Andrew Alliance is used with Pipette+.	OK <input type="checkbox"/> NOK <input type="checkbox"/>	Comment (if none, check this checkbox <input type="checkbox"/>)
Equip the power supply with the correct power adaptor based on the wall socket present in the lab.	OK <input type="checkbox"/> NOK <input type="checkbox"/>	Comment (if none, check this checkbox <input type="checkbox"/>)
Connect the power cord of the power supply to a compatible socket. Afterwards, plug the cable of the power supply into the back of Pipette+. 	OK <input type="checkbox"/> NOK <input type="checkbox"/>	Comment (if none, check this checkbox <input type="checkbox"/>)

<p>Verify that the “Plus” LED indication present on the front body of Pipette+ is blinking white. After few minutes, its color will turn blue or green.</p>	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>
---	---	--

D) Placement of pipettes on Pipette+ rack

<p>Verification that the right type of pipettes is installed on Pipette+</p>		
<p>Verify that the pipettes installed on Pipette+ stand are Andrew Alliance pipettes.</p>	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>
<p>Verify that the pipette display is oriented towards the user, as shown below.</p> 	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>
<p>Verify that the display color of each pipette placed in the rack of Pipette+, matches the color of the main button of the pipette (ex. Orange, blue, purple etc.) <u>Note:</u> this indicates that the pipette is properly connected and charging.</p>	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>

E) Pairing of Pipette+ to OneLab

<p>Association of Pipette+ to OneLab. This procedure is necessary only if Pipette+ hasn't been already paired to OneLab.</p>		
<p>If the Cloud version of OneLab is used, verify that the network, which Pipette+ will be connected to has Internet access and the webpage https://onelab.andrewalliance.com can be reached.</p>	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/> N/A <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>
<p>If the Cloud configuration of OneLab is used, make sure you have an OneLab account for the admin user that this device will be associated with. Visit https://onelab.andrewalliance.com to create a new free account.</p>	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>



<p>In case of a Standalone local server verify that the server hosting OneLab is placed preferably on the right side of Pipette+ and at least 3 meters away from the Pipette+ device (if used on Wi-Fi) or as far as the Ethernet cable allows in case of wired connection. <u>Note:</u> placing the server closer could create interference with the pipettes placed on the Pipette+.</p>	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/> N/A <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>
<p>In case of local server (Standalone) or Enterprise configuration, verify that the server hosting OneLab is turned ON and the network is available.</p>	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/> N/A <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>
<p>Follow the procedure of “Getting started” as explained in Pipette+ Quick Start Guide. After the pairing procedure is successful, Pipette+ should be visible in the Device page of OneLab having the status “Ready”.</p>	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>
<p>If Ethernet connection will be used with Pipette+ instead of Wi-fi, make sure you use the Ethernet cable provided by Andrew Alliance. <u>Note:</u> if a longer cable is required, the customer IT who provided the cable must make sure the cable is properly functional.</p> 	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/> N/A <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>
<p>In case of pairing to an OneLab Standalone or Enterprise server, make sure to specify the correct OneLab address during the pairing process:</p> <p><u>Standalone:</u> OneLab web address <input type="text" value="http://192.168.5.1"/></p> <p><u>Enterprise:</u> OneLab web address <input type="text" value="The OneLab Enterprise address provided by the IT"/></p> <p><u>Note:</u> the Standalone Server is equipped with two Ethernet ports: the one on top must be used in case of Ethernet connection.</p>	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/> N/A <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>

After the pairing is successfully completed, access your OneLab. In that moment, select the Lab you would like to link the Pipette+.	OK <input type="checkbox"/> NOK <input type="checkbox"/>	<u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)
Select the current Pipette+ in the list of available devices present in OneLab Device Tab. Click on “Locate Me” to verify that the communication between OneLab and Pipette+ is successful. Pipette+ should make a sound and its “Plus” LED should blink violet until you click again on “Locate Me”.	OK <input type="checkbox"/> NOK <input type="checkbox"/>	<u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)
Click on the button “Refresh tools” and verify that the pipettes present on Pipette+ stand are visible in OneLab.	OK <input type="checkbox"/> NOK <input type="checkbox"/>	<u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)

1.2.IQ test report

Below Installation Qualification Report is shown. The Installation Qualification is successful if the instrument passed **ALL** the tests. The report should be signed by both the person who executed the IQ process, as well as the owner of the device.

INSTALLATION QUALIFICATION REPORT				
<i>IQ test</i>	<i>DESCRIPTION</i>	<i>RESULT</i>	<i>COMMENT (if none, bar the cell)</i>	<i>Evaluation (Pass/Fail)</i>
A	User Manual and Documents	OK <input type="checkbox"/> NOK <input type="checkbox"/>		PASS <input type="checkbox"/> FAIL <input type="checkbox"/>
B	Software installation	OK <input type="checkbox"/> NOK <input type="checkbox"/>		PASS <input type="checkbox"/> FAIL <input type="checkbox"/>
C	Installation of Pipette+	OK <input type="checkbox"/> NOK <input type="checkbox"/>		PASS <input type="checkbox"/> FAIL <input type="checkbox"/>
D	Placement of pipettes on Pipette+	OK <input type="checkbox"/> NOK <input type="checkbox"/>		PASS <input type="checkbox"/> FAIL <input type="checkbox"/>
E	Power On the device	OK <input type="checkbox"/> NOK <input type="checkbox"/>		PASS <input type="checkbox"/> FAIL <input type="checkbox"/>
F	Pairing of Pipette+ to OneLab	OK <input type="checkbox"/> NOK <input type="checkbox"/>		PASS <input type="checkbox"/> FAIL <input type="checkbox"/>

	Function	Name	Date	Signature
Performed by				
Approved by				

2. Operational Qualification (OQ)

Document #: 618.002		Revision: 7	
	Function	Name	Date
Author	Lab Automation Service & Support Engineer	L.M. Turrado	25 Feb 2022
Reviewed by	Global Service & Support Sr Manager	N. Vamvoukas	25 Feb 2022

The **Operational Qualification (OQ)** is a test process that evaluates the proper functioning of a facility or an appliance. The Operational Qualification (OQ) may only be performed after a successfully completed Installation Qualification (IQ). The documentation for an Operational Qualification consists of:

- OQ test plan
- OQ report

2.1. OQ test plan

The **OQ test plan** contains detailed information of the tests to be performed on the respective equipment in the context of the Operation Qualification (OQ). This test plan contains detailed specifications on the course of the test itself.

After completion of the OQ test plan, the results are summarized, evaluated and combined in an OQ report. All deviations, as well as measures taken to eliminate these, are documented. After the deviations have been eliminated, they must be tested again, documented, and evaluated in a new OQ report.

The Operational Qualification template for Pipette+ is reported in the next page.

OPERATION QUALIFICATION			
System name:		Pipette+	
Serial number:		618. _____	
Document #	618.002	Revision	7

2.2. Getting started

This Operational Qualification procedure is based on the execution of a set of protocols to test the correct functionality of different modules of the system. Installation Qualification must be successfully performed before starting the Operation Qualification procedure.


It is suggested to use distilled water and a solution of Ponceau S with a concentration of 0.5 g/l. The use of dye is suggested to have a better visibility of the correctly filled wells and consumables.

As Pipette+ handles electronic pipettes, there is no need to test the fluidic performances of the system at the time of the installation. The latter have been tested during the manufacturing process of the pipettes. A performance qualification document is delivered along with the pipettes.

A) Experiment execution and verification correct execution

Select one of the options below depending on your Pipette setup.

OPTION 1		SELECTED <input type="checkbox"/>
Starter Kit configuration		N/A <input type="checkbox"/>
Verify that the following Andrew Alliance pipettes are installed on Pipette+ stand: <ul style="list-style-type: none"> - Single channel AA pipette 100-5000 µL - Single channel AA pipette 10-300 µL - 8-channel AA pipette 50-1200 µL 	OK <input type="checkbox"/> NOK <input type="checkbox"/>	<u>Comment (if none, check this checkbox <input type="checkbox"/>)</u>
Verify the availability of the following Sartorius pipette tips (one per type) for the execution of the experiment: <ul style="list-style-type: none"> - <u>300uL</u>: 790350, 790351, 790352, 790353, LH-L790350, LH-L790352, 790301F or LH-LF790301 - <u>1200uL</u>: 791210, 791211, 791212, 791213, LH-L791210, 791211F or LH-LF7991211 - <u>5mL</u>: 780304, 780305, 780300, 780308 	OK <input type="checkbox"/> NOK <input type="checkbox"/>	<u>Comment (if none, check this checkbox <input type="checkbox"/>)</u>
Verify the availability of the following labware for the execution of the experiment: <ul style="list-style-type: none"> - 50mL conical tube Corning #352070 or comparable to it. - Flat bottom 96-well plate Greiner™ #655101 or comparable to it. - 8-channel reservoir Integra #4332. - 1.5mL centrifuge conical microtube Eppendorf Safe-Lock #0030120086 or comparable to it. 	OK <input type="checkbox"/> NOK <input type="checkbox"/>	<u>Comment (if none, check this checkbox <input type="checkbox"/>)</u>
Execute the protocol "OQ Pipette+ v1.0" with Pipette+ and verify that all the pipettes in the list below have been selected (as shown on the picture below): <ul style="list-style-type: none"> - Single channel AA pipette 100-5000 µL - Single channel AA pipette 10-300 µL 	OK <input type="checkbox"/> NOK <input type="checkbox"/>	<u>Comment (if none, check this checkbox <input type="checkbox"/>)</u>

<p>- 8-channel AA pipette 50-1200 µL</p> 		
<p>During the experiment execution observe that in step #1 of the protocol:</p> <ul style="list-style-type: none"> • 100-5000 µL pipette has been selected. • A mixing step at 5mL volume has been done five times 	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/></p>	<p><i>Comment (if none, check this checkbox <input type="checkbox"/>)</i></p>
<p>During the experiment execution observe that in step #2 of the protocol:</p> <ul style="list-style-type: none"> • The same pipette tip pre-wetted in the previous step has been used • Repetitive dispensing of 1mL aliquot has been done in five different microtube 	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/></p>	<p><i>Comment (if none, check this checkbox <input type="checkbox"/>)</i></p>
<p>During the experiment execution observe that in step #3 of the protocol:</p> <ul style="list-style-type: none"> • A different pipette tip is used (tip ejection) • A volume of 5mL has been dispensed using reverse pipetting technique. • The excess volume has been purged back to the source 	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/></p>	<p><i>Comment (if none, check this checkbox <input type="checkbox"/>)</i></p>
<p>During the experiment execution observe that in step #4 of the protocol:</p> <ul style="list-style-type: none"> • 8-channel 50-1200 µL has been selected • A volume of 100 µL has been dispensed using forward pipetting technique from the reservoir to column 1 and 2 of the 96-well plate using an air cushion after the aspiration • For each dispensing, new pipette tips have been used 	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/></p>	<p><i>Comment (if none, check this checkbox <input type="checkbox"/>)</i></p>
<p>During the experiment execution observe that in step #5 of the protocol:</p> <ul style="list-style-type: none"> • 1-channel 10-300 µL has been selected • Mixing of 70µL in A1 well of the 96-well plate two times • Aspiration of 50µL volume from A1 of the 96-well plate • Dispensing of 50µL volume in A2 of the 96-well plate • Mixing of 100µL volume two times 	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/></p>	<p><i>Comment (if none, check this checkbox <input type="checkbox"/>)</i></p>

<p>During the experiment execution observe that in step #6 of the protocol:</p> <ul style="list-style-type: none"> • A different pipette tip has been used • 1-channel 10-300 µL has been selected • Mixing of 70µL at the in B1 well of the 96-well plate two times • Aspiration of 50µL volume from B1 of the 96-well plate • Dispensing of 50µL volume in B2 of the 96-well plate • Mixing of 100µL volume two times 	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>
<p>The protocol was executed without errors.</p>	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>

<p>OPTION 2 Custom Pipette configuration</p>		<p>SELECTED <input type="checkbox"/> N/A <input type="checkbox"/></p>
<p>Select one of the below configurations:</p> <ul style="list-style-type: none"> • Option 1: The pipettes present in the Pipette+ can handle volumes from 0.5uL up to 1200uL with single and/or multi-channel pipettes. • Option 2: The pipettes present in the Pipette+ can handle volumes from 10uL to 5mL with single ONLY or single/ multi-channel pipettes. • Option 3: The pipettes present in the Pipette+ can handle volumes from 10uL to 1200uL with multi-channel ONLY pipettes. • Option 4: The pipettes present in the Pipette+ can handle big volumes only (5mL and/or 10mL pipettes) 	<p>Option 1 <input type="checkbox"/> Option 2 <input type="checkbox"/> Option 3 <input type="checkbox"/> Option 4 <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>
<p>Execute the relevant operation test protocol depending on the option selected in the previous step:</p> <ul style="list-style-type: none"> • Option 1: "OQ custom from 0.5uL up to 1200uL Single and_or Multi.onp" • Option 2: "OQ custom 10uL to 5mL with single ONLY or Single Multichannel.onp " • Option 3: "OQ custom from 10uL to 1200uL multi-channel ONLY.onp " • Option 4: "OQ custom big volumes only 5mL and_or 10mL pipettes. onp" <p><u>Note:</u> It might be required to replace the labware used in the OQ test protocol with a labware of similar properties that you have available.</p>	<p>OK <input type="checkbox"/> NOK <input type="checkbox"/></p>	<p><u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)</p>

Follow the instructions stated in OneLab and perform the steps of the chosen protocol.	OK <input type="checkbox"/> NOK <input type="checkbox"/>	<u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)
The protocol was executed without errors.	OK <input type="checkbox"/> NOK <input type="checkbox"/>	<u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)

B) Collection of Logfiles

<p>Following the completion of the OQ protocol(s), collect the logfiles.</p> <p>Important note: the collection of logfiles must be performed right after the protocols have been completed and without power cycling the device. If for any reason the unit had to be power cycled, the OQ process and protocol execution must be restarted from the beginning.</p>	OK <input type="checkbox"/> NOK <input type="checkbox"/>	<u>Comment</u> (if none, check this checkbox <input type="checkbox"/>)
---	---	---

2.3.OQ test report

Below Operation Qualification Report is shown. The Operation Qualification is successful if the instrument passed **ALL** the tests. The report must be signed by both the person who executed the OQ process, as well as the owner of the device.

OPERATION QUALIFICATION REPORT				
OQ test	DESCRIPTION	RESULT	COMMENT (if none, bar the cell)	Evaluation (Pass/Fail)
A	Experiment execution and verification of correct operation	OK <input type="checkbox"/> NOK <input type="checkbox"/>		PASS <input type="checkbox"/> FAIL <input type="checkbox"/>
B	Collection of logfiles	OK <input type="checkbox"/> NOK <input type="checkbox"/>		PASS <input type="checkbox"/> FAIL <input type="checkbox"/>

	Function	Name	Date	Signature
Performed by				
Approved by				

3. Contact

For any questions or information, do not hesitate to:

- Visit our website www.andrewalliance.com
- Contact us through our Intercom messenger present in [OneLab](#)
- Send us an email at aa_support@waters.com